


**Alligator™**  
Retrieval Device

# Navigate. Activate. Eliminate.

Get A Grip On More  
Confident Retrievals.



**ALLIGATOR™**—The industry's leading retrieval device for neurovascular procedures.

## Putting Confident Control And End-To-End Performance Close At Hand.

The Alligator™ Retrieval Device provides physicians with an advanced solution, specifically designed and approved for neurovascular procedures.

### Proximal Access, For Procedural Confidence

By allowing access from the front side of the foreign body—rather than requiring the additional time and perforation risks involved with going distal—Alligator is designed to optimize safety, accelerate removal, minimize risk, and maximize successful outcomes.

### Knows When To Hold, Knows When To Fold

The four interlocking jaws deliver powerful and consistent gripping strength at a distance, while our “catch and release” design provides unsurpassed control, the ability to reposition foreign bodies for efficient removal, and instant adaptation to changing procedural needs.

### Micro-Movement Precision


Constructed of radiopaque platinum tip and a stainless steel core wire, the Alligator delivers clear visualization under fluoroscopy, remarkably tactile feedback, “as-expected” movement control, and unsurpassed precision through even the most intricate anatomies.

### Atraumatic Design And Ease Of Use

Alligator's atraumatic gripping arms are easily deployed by simply advancing or retreating the microcatheter when adjacent to the target—WITHOUT additional levers or controls. This allows for positive and direct grasping, engagement of the foreign body at any point, “hands-on” sensitivity, and optimized retrieval during critical procedures.

### Versatile Compatibility

Engineered for superb guidewire flexibility and trackability, Alligator's range of sizes allows for precise parent vessel matching, and maximum vessel engagement. And it can be used with any 3F (.021") microcatheter—making for easy compatibility.



### The Easy Way To Get Out Of Tough Situations.


- Proximal access for speed and safety
- Sensitive, tactile feedback
- Precise micro-movement control
- Solid gripping strength at a distance
- “Catch and release” ability for instant repositioning

## Alligator Retrieval Device

Catalog Number	Description	Quantity
FA-88810-20	2 mm ARD	1
FA-88810-30	3 mm ARD	1
FA-88810-40	4 mm ARD	1
FA-88810-50	5 mm ARD	1

### Multiple sizes to optimize retrieval

Alligator Retrieval Device (ARD) sized to parent vessel completely opens for maximum engagement.



### Indication for use

This device is intended for use in the peripheral and neurovasculature for the retrieval of foreign objects.

### Contra Indications:

This device is not intended for:

- Removal of foreign objects entrapped by tissue growth
- Removal of implanted pacing leads

### Potential Complications

Potential complications include, but are not limited to, the following:

- Vessel perforation or dissection
- Thromboemboli
- Hemorrhage
- Vasospasm
- Neurological deficits including stroke and death

See package insert for complete indications, contra indications, potential complications, warnings, and instructions for use.

#### NEUROVASCULAR | PERIPHERAL VASCULAR

Access · Balloons · Carotid · Embolic Coils · Embolic Protection · Flow Diversion · Flow Restoration · Liquid Embolics · Plaque Excision · Procedural Support  
PTA Balloons · Remodeling Devices · Retrieval Devices · Stents · Thrombectomy

ev3 Europe  
International Headquarters  
106-108 rue La Boétie  
75008 Paris  
France  
PH +33 156 88 59 10  
FX +33 156 88 59 11

ev3 Corporate  
World Headquarters  
Peripheral Vascular  
9600 54th Avenue North  
Plymouth, MN 55442-2111  
USA  
PH +1 763 398 7000  
FX +1 763 398 7001  
www.ev3.net

ev3 Neurovascular  
9775 Toledo Way  
Irvine, CA 92618 USA  
PH +1 949 837 3700  
FX +1 949 837 2044  
  
ev3 International  
Distribution Centre  
Europalaan 25  
6199 AB Maastricht-Airport  
The Netherlands  
PH +31 (0) 433 659 220  
FX +31 (0) 43 364 6395

ev3 SAS France  
PH +33 (0) 156 88 31 10  
FX +33 (0) 156 88 31 11

ev3 B.V. Benelux  
PH +31 (0) 433 659 223  
FX +31 (0) 433 650 283

ev3 Technologies  
Iberica, S.L. Spain  
PH +34 91 656 7154  
FX +34 91 656 7214

ev3 S.r.l. Italy  
PH +39 0267 977 61  
FX +39 0266 711 637

ev3 Nordic AB  
PH +46 859 000 950  
FX +46 859 000 959

ev3 GmbH Germany, Austria  
PH +49 228 528 830  
FX +49 228 528 8360

ev3 Ltd. United Kingdom  
PH +44 1279 659 900  
FX +44 1279 654 900

