

NanoCross™

.014 PTA Balloon

Now Available in

210 mm
length

Reach your
destination.

The next generation
of .014 PTA Balloons



Introducing NanoCross™




The next generation in PTA Balloon technology from ev3.

Reach your destination with an enhanced OTW

.014 balloon designed to cross smoothly and deflate quickly.

Ultra-low Crossing Profile


- ThermoFlex™ Balloon Material designed to optimize thin-wall strength
- SlimTec™ Balloon Folding Process designed to provide lowest .014 crossing profile



3 x 150 NanoCross, 3 x 120 Competitor A and 3 x 120 Competitor B

Improved Deliverability

- Robust catheter design improves pushability and tracking
- CrossCoat™ Hydrophilic Balloon Coating reduces friction



* ev3 data on file




Now Available in 210 mm Tapered Length

- The balloons taper .5 mm from proximal to distal

Marker Bands

- Two additional marker bands centered on the 210 mm balloon provide additional visibility


Greatly Reduced Deflation Times



Bevel 360™

Revolutionary Bevel 360™ Tip Design

- 360 degree bevel provides ultra smooth transition from wire to tip
- Lowest tip entry profile available



NanoCross™ .014 PTA Balloon

Product Code 150 cm catheter length	Product Code 90 cm catheter length	Balloon Diameter (mm)	Balloon Length (mm)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Recommended Introducer Sheath (F)
AB14W015020150	AB14W015020090	1.5	20	10	14	4
AB14W020020150	AB14W020020090	2.0	20	10	14	4
AB14W020040150	AB14W020040090	2.0	40	10	14	4
AB14W020080150	AB14W020080090	2.0	80	10	14	4
AB14W020120150	AB14W020120090	2.0	120	10	14	4
AB14W020150150	AB14W020150090	2.0	150	10	14	4
AB14W020210150	AB14W020210090	2.0 (proximal) / 1.5 (distal)	210	10	14	4
AB14W025020150	AB14W025020090	2.5	20	10	14	4
AB14W025040150	AB14W025040090	2.5	40	10	14	4
AB14W025080150	AB14W025080090	2.5	80	10	14	4
AB14W025120150	AB14W025120090	2.5	120	10	14	4
AB14W025150150	AB14W025150090	2.5	150	10	14	4
AB14W025210150	AB14W025210090	2.5 (proximal) / 2.0 (distal)	210	10	14	4
AB14W030020150	AB14W030020090	3.0	20	10	14	4
AB14W030040150	AB14W030040090	3.0	40	10	14	4
AB14W030080150	AB14W030080090	3.0	80	10	14	4
AB14W030120150	AB14W030120090	3.0	120	10	14	4
AB14W030150150	AB14W030150090	3.0	150	10	14	4
AB14W030210150	AB14W030210090	3.0 (proximal) / 2.5 (distal)	210	10	14	4
AB14W035020150	AB14W035020090	3.5	20	8	14	4
AB14W035040150	AB14W035040090	3.5	40	8	14	4
AB14W035080150	AB14W035080090	3.5	80	8	14	4
AB14W035120150	AB14W035120090	3.5	120	8	14	4
AB14W035150150	AB14W035150090	3.5	150	8	14	4
AB14W035210150	AB14W035210090	3.5 (proximal) / 3.0 (distal)	210	8	14	4
AB14W040020150	AB14W040020090	4.0	20	7	14	4
AB14W040040150	AB14W040040090	4.0	40	7	14	4
AB14W040080150	AB14W040080090	4.0	80	7	14	4
AB14W040120150	AB14W040120090	4.0	120	7	14	4
AB14W040150150	AB14W040150090	4.0	150	7	14	4
AB14W040210150	AB14W040210090	4.0 (proximal) / 3.5 (distal)	210	8	14	4

Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

NanoCross, ThermoFlex, SlimTec and CrossCoat are all trademarks of ev3 Inc. Other names appearing in this document are the property of their respective owners.

NEUROVASCULAR | PERIPHERAL VASCULAR

Access · Balloons · Embolic Coils · Embolic Protection · Flow Diversion · Flow Restoration · Liquid Embolics · Plaque Excision · Procedural Support · Remodeling Devices · Retrieval Devices · Stents

ev3 Europe
International Headquarters
106-108 rue La Boétie
75008 Paris
France
PH +33 156 88 59 10
FX +33 156 88 59 11

ev3 Corporate
World Headquarters
Peripheral Vascular
3033 Campus Drive
Plymouth, MN 55441
USA
PH +1 763 398 7000
FX +1 763 398 7001
TOOL-FREE +1 800 716 6700

ev3 Neurovascular
9775 Toledo Way
Irvine, CA 92618 USA
PH +1 949 837 3700
FX +1 949 837 2044

ev3 International
Distribution Centre
Europalaan 25
6199 AB Maastricht-Airport
The Netherlands
PH +31 (0) 433 659 220
FX +31 (0) 43 364 6395

ev3 SAS France
PH +33 (0) 156 88 31 10
FX +33 (0) 156 88 31 11

ev3 B.V. Benelux
PH +31 (0) 433 659 223
FX +31 (0) 433 650 283

ev3 Technologies
Iberica, S.L. Spain
PH +34 91 656 7154
FX +34 91 656 7214

ev3 S.r.l Italy
PH +39 0267 977 61
FX +39 0266 711 637

ev3 Nordic AB
PH +46 859 000 950
FX +46 859 000 959

ev3 Sp z o.o. Poland
PH +48 32 747 01 44
FX +48 32 747 01 45

ev3 GmbH Germany, Austria
PH +49 228 528 830
FX +49 228 528 8360

ev3 Ltd. United Kingdom
PH +44 1279 659 900
FX +44 1279 654 900

www.ev3.net

