

## Stents

### Protégé™ GPS™

Self-Expanding Peripheral  
Stent System

# The ILIAC *solution*

### Visi-Pro™


Balloon-Expandable  
Peripheral Stent System

## Visi-Pro™

0.035" over the wire system: *Balloon-Expandable Peripheral Stent System*



Effortless  
*precision*



## Protégé™ GPS™

0.035" over the wire system: *Self-Expanding Peripheral Stent System*

Superb  
*control*

## Accurate Placement

- Minimal stent shortening
- Radiopaque markers at both ends of the stent
- Ideal for ostial placement
- Minimal balloon to stent overhang

## Maximum Strength

- Excellent radial strength and support
- Optimal cell design for conformability and trackability

## Superior Deliverability and Trackability

- Low profile
- Broad offering of 6 Fr



---

## Meeting Every Requirement for Iliac Stenting

- Optimal expansion force
- High flexibility
- Full line of diameters from 6 to 14 mm

## Ensuring Superb Control and Accurate Deployment

- Tantalum GPS Markers on ends of stent
- Unique EX.P.R.T.™ system eliminates jumping
- Minimal stent shortening

## Achieving Fast Access

- Atraumatic tip
- Flexible catheter
- All sizes 6 Fr sheath compatible



## Visi-Pro™

### 0.035" Balloon-Expandable Peripheral Stent System

#### ORDERING INFORMATION

Product Number Catheter Length 80 cm	Product Number Catheter Length 135 cm	Stent Diameter (mm)	Stent Length (mm)	Balloon Length (mm)	Sheath Size (Fr)
PXP35-05-12-080	-	5	12	15	6*
PXP35-05-17-080	PXP35-05-17-135	5	17	20	6
PXP35-05-27-080	PXP35-05-27-135	5	27	30	6
PXP35-05-37-080	PXP35-05-37-135	5	37	40	6
PXP35-05-57-080	PXP35-05-57-135	5	57	60	6
PXP35-06-12-080	-	6	12	15	6
PXP35-06-17-080	PXP35-06-17-135	6	17	20	6
PXP35-06-27-080	PXP35-06-27-135	6	27	30	6
PXP35-06-37-080	PXP35-06-37-135	6	37	40	6
PXP35-06-57-080	PXP35-06-57-135	6	57	60	6
PXP35-07-12-080	-	7	12	15	6
PXP35-07-17-080	PXP35-07-17-135	7	17	20	6
PXP35-07-27-080	PXP35-07-27-135	7	27	30	6
PXP35-07-37-080	PXP35-07-37-135	7	37	40	6
PXP35-07-57-080	PXP35-07-57-135	7	57	60	6
PXP35-08-17-080	PXP35-08-17-135	8	17	20	6
PXP35-08-27-080	PXP35-08-27-135	8	27	30	6
PXP35-08-37-080	PXP35-08-37-135	8	37	40	6
PXP35-08-57-080	PXP35-08-57-135	8	57	60	6
PXP35-09-17-080	PXP35-09-17-135	9	17	20	7
PXP35-09-27-080	PXP35-09-27-135	9	27	30	7
PXP35-09-37-080	PXP35-09-37-135	9	37	40	7
PXP35-09-57-080	PXP35-09-57-135	9	57	60	7
PXP35-10-17-080	PXP35-10-17-135	10	17	20	7
PXP35-10-27-080	PXP35-10-27-135	10	27	30	7
PXP35-10-37-080	PXP35-10-37-135	10	37	40	7
PXP35-10-57-080	PXP35-10-57-135	10	57	60	7

\*6 Fr=0.085" I.D.

Visi-Pro carries the CE mark for use in the renal, common and external iliac, and subclavian arteries, as well as malignant biliary use. Visi-Pro is a trademark of ev3 Inc.

Protégé GPS carries the CE mark for use in the common iliac, external iliac, subclavian arteries, and malignant biliary use. Protégé and GPS are trademarks of ev3 Inc.

Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

## Protégé™ GPS™

### 0.035" Self-Expanding Peripheral Stent System

#### ORDERING INFORMATION

Product Number Catheter Length 80 cm	Product Number Catheter Length 120 cm	Stent Diameter (mm)	Stent Length (mm)	Lumen Diameter (mm)	Sheath Size (Fr)
SERP65-06-20-80	SERP65-06-20-120	6	20	4.5-5.5	6
SERP65-06-30-80	SERP65-06-30-120	6	30	4.5-5.5	6
SERP65-06-40-80	SERP65-06-40-120	6	40	4.5-5.5	6
SERP65-06-60-80	SERP65-06-60-120	6	60	4.5-5.5	6
SERP65-06-80-80	SERP65-06-80-120	6	80	4.5-5.5	6
SERP65-07-20-80	SERP65-07-20-120	7	20	5.5-6.5	6
SERP65-07-30-80	SERP65-07-30-120	7	30	5.5-6.5	6
SERP65-07-40-80	SERP65-07-40-120	7	40	5.5-6.5	6
SERP65-07-60-80	SERP65-07-60-120	7	60	5.5-6.5	6
SERP65-07-80-80	SERP65-07-80-120	7	80	5.5-6.5	6
SERP65-08-20-80	SERP65-08-20-120	8	20	6.5-7.5	6
SERP65-08-30-80	SERP65-08-30-120	8	30	6.5-7.5	6
SERP65-08-40-80	SERP65-08-40-120	8	40	6.5-7.5	6
SERP65-08-60-80	SERP65-08-60-120	8	60	6.5-7.5	6
SERP65-08-80-80	SERP65-08-80-120	8	80	6.5-7.5	6
SERP65-09-20-80	SERP65-09-20-120	9	20	7.5-8.5	6
SERP65-09-30-80	SERP65-09-30-120	9	30	7.5-8.5	6
SERP65-09-40-80	SERP65-09-40-120	9	40	7.5-8.5	6
SERP65-09-60-80	SERP65-09-60-120	9	60	7.5-8.5	6
SERP65-09-80-80	SERP65-09-80-120	9	80	7.5-8.5	6
SERP65-10-20-80	SERP65-10-20-120	10	20	8.5-9.5	6
SERP65-10-30-80	SERP65-10-30-120	10	30	8.5-9.5	6
SERP65-10-40-80	SERP65-10-40-120	10	40	8.5-9.5	6
SERP65-10-60-80	SERP65-10-60-120	10	60	8.5-9.5	6
SERP65-10-80-80	SERP65-10-80-120	10	80	8.5-9.5	6
SERP65-12-20-80	SERP65-12-20-120	12	20	9.5-11.0	6
SERP65-12-30-80	SERP65-12-30-120	12	30	9.5-11.0	6
SERP65-12-40-80	SERP65-12-40-120	12	40	9.5-11.0	6
SERP65-12-60-80	SERP65-12-60-120	12	60	9.5-11.0	6
SERP65-12-80-80	SERP65-12-80-120	12	80	9.5-11.0	6
SERP65-14-20-80	SERP65-14-20-120	14	20	11.5-13.0	6
SERP65-14-30-80	SERP65-14-30-120	14	30	11.5-13.0	6
SERP65-14-40-80	SERP65-14-40-120	14	40	11.5-13.0	6
SERP65-14-60-80	SERP65-14-60-120	14	60	11.5-13.0	6
SERP65-14-80-80	SERP65-14-80-120	14	80	11.5-13.0	6

#### VISI-PRO COMPLIANCE CHART

Visi-Pro Diameter (mm)	Inflation Pressure (atm)				
	8	9	10	11	12
5.0	5.00	5.09	5.16	5.22	5.28
6.0	6.00	6.11	6.22	6.31	6.39
7.0			7.00	7.09	7.17
8.0			8.00	8.15	8.26
9.0			9.00	9.15	9.28
10.0			10.00	10.11	10.21

● Diameter at Nominal Pressure ● Diameter at Rated Burst Pressure

#### NEUROVASCULAR | PERIPHERAL VASCULAR

Access · Balloons · Embolic Coils · Embolic Protection · Flow Diversion · Flow Restoration · Liquid Embolics · Plaque Excision · Procedural Support · Remodeling Devices · Retrieval Devices · Stents

**ev3 Europe**  
International Headquarters  
106-108 rue La Boétie  
75008 Paris  
France  
PH +33 156 88 59 10  
FX +33 156 88 59 11

**ev3 Corporate**  
World Headquarters  
Peripheral Vascular  
3033 Campus Drive  
Plymouth, MN 55441  
USA  
PH +1 763 398 7000  
FX +1 763 398 7001

**ev3 Neurovascular**  
9775 Toledo Way  
Irvine, CA 92618 USA  
PH +1 949 837 3700  
FX +1 949 837 2044

**ev3 International**  
Distribution Centre  
Europalaan 25  
6199 AB Maastricht-Airport  
The Netherlands  
PH +31 (0) 433 659 220  
FX +31 (0) 43 364 6395

**ev3 SAS France**  
PH +33 (0) 156 88 31 10  
FX +33 (0) 156 88 31 11

**ev3 B.V. Benelux**  
PH +31 (0) 433 659 223  
FX +31 (0) 433 650 283

**ev3 Technologies**  
Iberica, S.L. Spain  
PH +34 91 656 7154  
FX +34 91 656 7214

**ev3 S.r.l Italy**  
PH +39 0267 977 61  
FX +39 0266 711 637

**ev3 Nordic AB**  
PH +46 859 000 950  
FX +46 859 000 959

**ev3 Sp z o.o. Poland**  
PH +48 32 747 01 44  
FX +48 32 747 01 45

**ev3 GmbH Germany, Austria**  
PH +49 228 528 830  
FX +49 228 528 8360

**ev3 Ltd. United Kingdom**  
PH +44 1279 659 900  
FX +44 1279 654 900

www.ev3.net

