


**TrailBlazer™**




Procedural Support  
Catheter

Crossing Boundaries.  
Expanding Options.




# TrailBlazer™ Procedural Support Catheter

Guidewire Compatibility	Marker Band Spacing (mm)	Working Length			
		65 cm	90 cm	135 cm	150 cm
0.014"	15			SC-014-135	SC-014-150
0.018"	15		SC-018-090	SC-018-135	SC-018-150
0.035"	50	SC-035-065	SC-035-090	SC-035-0135	SC-035-150



## Increased Pushability

- Unique PushPlus™ Catheter design provides an average of 29% more pushability than the leading competitive support catheters
- Robust outer layer catheter provides superior pushability
- Lubricious inner layer catheter enables smooth guidewire movement



## Tapered, Low Profile Tip for Ease of Lesion Entry

- Best in class lesion entry profile on 0.014" and 0.018"
- Rounded tip provides ultra smooth transition to wire
- Increased pushability without sacrificing profile

## Visible, Seamless Catheter

- Three embedded marker bands provide seamless catheter; essentially eliminates interaction with previously implanted devices
- CrossCoat™ Hydrophilic Coating on distal 40 cm for smooth tracking through challenging vasculature
- 50% longer marker bands designed for enhanced visibility; 1.5 mm vs. competitive 1 mm

Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. TrailBlazer™ is trademark of ev3 Inc.

### NEUROVASCULAR | PERIPHERAL VASCULAR

Access · Balloons · Embolic Coils · Embolic Protection · Flow Diversion · Flow Restoration · Liquid Embolics · Plaque Excision · Procedural Support · Remodeling Devices · Retrieval Devices · Stents

**ev3 Europe International Headquarters**  
106-108 rue La Boétie  
75008 Paris  
France  
PH +33 156 88 59 10  
FX +33 156 88 59 11

**ev3 Corporate World Headquarters Peripheral Vascular**  
3033 Campus Drive  
Plymouth, MN 55441  
USA  
PH +1 763 398 7000  
FX +1 763 398 7001

**ev3 Neurovascular**  
9775 Toledo Way  
Irvine, CA 92618 USA  
PH +1 949 837 3700  
FX +1 949 837 2044

**ev3 International Distribution Centre**  
Europalaan 25  
6199 AB Maastricht-Airport  
The Netherlands  
PH +31 (0) 433 659 220  
FX +31 (0) 43 364 6395

**ev3 SAS France**  
PH +33 (0) 156 88 31 10  
FX +33 (0) 156 88 31 11

**ev3 B.V. Benelux**  
PH +31 (0) 433 659 223  
FX +31 (0) 433 650 283

**ev3 Technologies Iberica, S.L. Spain**  
PH +34 91 656 7154  
FX +34 91 656 7214

**ev3 S.r.l. Italy**  
PH +39 0267 977 61  
FX +39 0266 711 637

**ev3 Nordic AB**  
PH +46 859 000 950  
FX +46 859 000 959

**ev3 Sp z o.o. Poland**  
PH +48 32 747 01 44  
FX +48 32 747 01 45

**ev3 GmbH Germany, Austria**  
PH +49 228 528 830  
FX +49 228 528 8360

**ev3 Ltd. United Kingdom**  
PH +44 1279 659 900  
FX +44 1279 654 900

www.ev3.net

EC REP  
MedMark Europe Sarl  
BP 2332  
F-38033 Grenoble Cedex 2  
FRANCE

